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EXEMPLO:

Legal duty to inform about battery disposal

Batteries do not belong in household waste. As a consumer you are legally obligated to return batteries. You can hand in your old batteries at the public collection point in your community or anywhere where batteries of the same kind are sold.

Please note:

These symbols can be found on batteries:

Li = Battery contains lithium

Al = Battery contains alkali

Mn = Battery contains manganese

CR (Li); A (Al, Mn); AA (Al, Mn); AAA (Al, Mn)

Disposal of used electrical and electronic devices

This symbol on the product or its packaging indicates that the product is not to be treated like normal household waste, but must be handed over at a collection point for the recycling of electrical and electronic devices. Further information can be obtained in your community, at communal disposal services or the business where you bought the product.

Yours faithfully,

ADE (GmbH & Co.)

✂-----

Warranty section – 3 years on the body analysis scales

• Sender

• Reason for Complaint

• Date of purchase

Please send this section incl. proof of purchase back to your specialist retailer.

BODY FAT

BODY WATER

MUSCLE MASS

BONE MASS

ADE

Model Bella / BA 1402 Operating Manual

Table of Contents

1. How do the body analysis scales work?	Page 2
2. How to use them.....	Page 2
3. Start-up	Page 3
4. Programming / Inputting personal data	Page 4
5. Weighing function (weight measuring only)	Page 4
6. Body analysis	Page 4
7. Error messages	Page 7
8. Technical data	Page 7
9. Warranty	Page 7

Thank you for deciding to purchase this set of high-quality body analysis scales. Please read through this operating manual carefully before starting-up the scales for the first time, and keep it in a safe place so that this information is available to you when you need it.

1. How do the body analysis scales work?

With these new, revolutionary scales you will obtain information about your body composition in terms of weight, proportion of body fat (%), proportion of body water (%), muscle mass (%) and bone mass (in kg). Your personal data, such as sex (male or female), age and height, must first be entered into one of the individual storage cells. Taking a measurement is then as easy as weighing yourself on a completely normal set of scales. You need only choose your personal storage cell, stand on the scales and within a few seconds your current body status will be displayed for you. For measuring, a weak current is used, which gauges the resistance. Water, which is principally stored in the musculature, conducts very well; fat, on the other hand, does not. With reference to the measured resistance and the stored data, the body analysis scales can calculate your body composition very accurately.

2. How to use

A measurement should preferably be taken while you are undressed and barefoot.

Since the body is naturally subject to variations (e.g. through dehydration during sport or after a visit to the sauna, the intake of food or fluids), it is important to take measurements under the most consistent conditions possible, in order to obtain comparable values.

Care should be taken to ensure that measurements are always taken at the same time of day. In the morning, immediately after getting up, body water has still not distributed itself throughout the body. This means that the resistance is relatively high because the conductivity in the body is still quite poor. With increasing movement throughout the day the water distributes itself throughout the body, the conductivity of the body increases, the resistance and thus the proportion of fat are registered less. In order to determine the actual change of body composition it is therefore important to measure yourself always at the same time. The easiest way to achieve this is to take your measurements in the morning, ¼ hour after getting up either always before or always after you go to the toilet.

Despite the minimal current strength, people with pacemakers or other implants that transmit signals must not use these body analysis scales! With implants such as a titanium plate in the knee there is no danger. Nevertheless, metals do conduct currents extremely well, as a consequence of which the body fat value is “embellished”, thus the actual amount of fat is higher than that displayed on the scales. Trends in body composition can nevertheless be determined equally well. In women, periodic influences affect the results of measurements because hormonal variations influence the water balance of the body. If there is a lot of water in the body the fat values come out lower; if

7. Error messages

- “EEEE” = Overload indicator – the maximum capacity of the scales (180 kg) has been exceeded. Remove the load from the scales immediately, otherwise the load cell could incur damage.
Or: Faulty weighing or measuring process – e.g. when the scales are not standing stably or you are not standing still on the weighing platform.
- “Lo” = Batteries are empty – please replace the batteries with the battery type intended for your scales (4 x 1.5V AAA).
- “Err” = The body fat proportion is either too low (below 5% cannot be shown) or the body fat proportion is too high (over 50% cannot be shown).

8. Technical Data

Weight-bearing capacity x division:	Max. 180 kg x 100 g
Indicator resolution:	Proportion of body fat: 0.1%
Indicator resolution:	Proportion of body water: 0.1%
Indicator resolution:	Proportion of muscle mass: 0.1 %
Indicator resolution:	Bone mass: 100 g
Personal storage cells:	10
Age input:	10 – 85 years
Input height:	75 – 225 cm
Dimensions:	approx. 305 mm x 305 mm x 22 mm
Batteries:	4 x 1.5V AAA (contained in product package)

We reserve the right to make technical modifications

9. Warranty

ADE guarantees the free correction of defects resulting from material or manufacturing faults through repair or exchange 3 years from the purchase date. Please get your warranty section filled out and stamped by the retailer when you make your purchase. In case of warranty, please return the scales to your retailer along with the warranty section indicating the reason for your complaint.

CE compliance. This device suppresses electromagnetic interference in accordance with the applicable EC directive 2004/108/EG

Please note: Under extreme electromagnetic influences e.g. when operating a radio device in immediate proximity to the scales, the indicator value may be affected. After the disturbing influence has stopped, the product may once again be used as intended, although a re-activation may be necessary.

The proportion of body fat gives the amount of fat as a percentage of the total weight of the body. A blanket recommendation for this value does not exist, since this depends significantly on sex (because of different physiques) and age. Body fat is vital for day-to-day body function. It protects organs, cushions joints, regulates body temperature, stores vitamins and serves the body as an energy store. For women, a body fat proportion between 10% and 13% is essential; for men between 2% and 5% with the exception of competitive sportsmen.

The proportion of body fat is not an unambiguous indicator of health. Excessive body weight and a too-high proportion of body fat are being linked with the appearance of many diseases of affluence such as diabetes, cardiovascular illnesses, etc. Physical inactivity and poor diet are often triggers for these diseases. There is a clear causal connection between excess weight and physical inactivity.

A considerable part of the human body consists of water. This body water is variably distributed throughout the whole body. Lean muscle tissue contains about 75% water, blood around 83% water, body fat approx. 25% water and bones have a water proportion of about 22%.

In men approx. 60% of body mass consists of water. In women, that value lies around 55% (determined by a higher proportion of body fat). Of this total body water approx. two thirds is stored in the cells and is therefore called intracellular water. The other third is extracellular water.

Women 					
Age	Body fat				Body water
	Slender	Normal	Slightly overweight	Overweight	Normal
10-16	<18%	18-28%	29-35%	>35%	57-67%
17-39	<20%	20-32%	33-38%	>38%	47-57%
40-55	<23%	23-35%	36-41%	>41%	42-52%
56-85	<24%	24-36%	37-42%	>42%	37-47%
Men 					
Age	Body fat				Body water
	Slender	Normal	Slightly overweight	Overweight	Normal
10-16	<10%	10-18%	19-23%	>23%	58-72%
17-39	<12%	12-20%	21-25%	>25%	53-67%
40-55	<13%	13-21%	22-26%	>26%	47-61%
56-85	<14%	14-22%	23-27%	>27%	42-56%

For muscle mass there are no universal standard values.

there is less water, then somewhat higher. For pregnant women there is no danger to their unborn children. Nevertheless, during pregnancy the values are so distorted by the water retention that the measurement has little meaning.

The ADE company assumes no liability for damages or losses caused by the body analysis scales, nor for third-party claims. This product is intended exclusively for private use. It is not intended for professional use in hospitals or medical institutions.

3. Start-up

Safety information:

This device is not intended for use by people (including children) with limited physical, sensory or mental capacities or lacking experience and knowledge, except when they are supervised by persons who are responsible for their safety or who are instructed in the use of the device. Children should not be allowed to play with the device.

Not suitable for people with electronic implants (pacemakers, etc.)

Place the scales on an even and firm surface (tiles, floorboards, etc.). Use on carpeted floors may result in incorrect measurements.

Attention! There is a danger of slipping when the surface of the scales is wet or when your feet are wet.

Please take care not to let any objects fall onto the scales, otherwise the glass could shatter.

Only clean the scales with a damp cloth. Do not use any solvents or abrasive cleaning products. Do not immerse in water.

Preparation:

Open the battery compartment cover on the underside of the scales and insert the provided batteries (4 x 1.5V AAA). Please be mindful of the correct polarity of the batteries (+/-).

On the underside of the scales, above the battery compartment, there is a small button for a country-specific unit setting (kg = kilogram / lb = pound / st = stone). Please make sure that the scales are set to the weight-unit relevant to you.

1. Programming / Inputting personal data

Prior to first use your personal data such as sex (male or female), age and height must be entered and stored in one of the storage cells P0 – P9 (for up to 10 people). The scales make use of highly sensitive sensor buttons (Sensor Touch). The touch keys are located below the display. To enter your personal data use the buttons: SET (power-on and memory button), ▼ (up button), ▲ (down button) and follow the individual points:

1. Gently touch the SET button. "P0" will blink on the display. To choose the desired storage cell (P0 – P9), touch the ▲- or ▼- button and then save by touching the SET button.
2. After you have chosen the storage cell (e.g. storage cell P3), the symbol for sex will blink. To choose the desired sex (♂ male or ♀ female), touch the ▲- or ▼- button and then save by touching the SET button.
3. After you have chosen your sex (e.g. ♀ female), the indicator for age will blink on the display (default age 30). In order to enter your age, touch the ▲- or ▼- button and then save by touching the SET button.
4. After you have entered your age, the indicator for height will blink on the display (default 170 cm). In order to enter your height, touch the ▲- or ▼-button and then save by touching the SET button.

After you have finished entering your data, "0.0 kg" will appear on the display. The scales are now ready and you can undertake a body analysis. Without further activity the scales deactivate automatically after a few seconds. Your personal data has been saved.

If you touch the SET button again within 2 seconds after the first input, you can immediately enter the data of other people, and for a new entry you will repeat completely the steps mentioned above.

2. Weighing function (weight measurement only)

1. The scales have a step-on function. Simply step onto the scales, keep still and do not hold onto anything.
2. After the scales have stabilised, your weight will be flashed at you three times in succession and then be shown to you steadily for approximately 10 seconds.
3. After weighing, step off the scales. After a few seconds the scales will deactivate automatically.

3. Body analysis

An analysis is only possible when you are barefoot. The personal data must be entered beforehand (see page 4).

1. Activate the scales by touching the SET button. On the display the personal data of the last setting will appear first. If this is not your data, choose the storage field under which you saved your personal data using the ▲- or ▼-button. Your personal data will be shown. Confirm each of these by touching the SET button, after which the scales will return to "0.0 kg".
2. Now stand carefully on the scales, barefoot with your feet left and right on the stainless steel strips (measuring units), and remain still.
3. First you will be shown your weight. After that, while the body analysis is being carried out, the signal "o" will run on the display. Remain standing still. After the measurement has been taken the measured values will be shown to you on the display:

Your body weight in kg

1. Your proportion of body fat in % (FAT) ***
2. Your body water in % (TBW)
3. Your muscle mass in % (MUS)
4. Your bone mass in kg (BONE)

*** While the proportion of body fat is being shown, at the very bottom of the display appears an evaluation:

	= fat proportion too low,
	= fat proportion in healthy range
	= fat proportion elevated
	= fat proportion too high

After that your measured values (weight, body fat, body water, muscle mass and bone mass) will be shown to you in turn, 3 times in succession.

EXAMPLE: