Körperanalysewaage Emmy

Body Analyser Scale Emmy | Báscula de análisis corporal Emmy |
Balance d'analyse corporelle Emmy | Bilancia con analisi corporea
Emmy | Waga z analizą organizmu Emmy | Lichaamsanalyseschaal
Emmy

BA 2500-1

BA 2500-2

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BA 2500-1 | BA 2500-2

Operating instructions | Instrucciones de operación | Mode d'emploi
Istruzioni per l'uso | Instrukcja obsługi | Gebruiksaanwijzing

Liebe Kundin, lieber Kunde !

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Dear Customer!

You have chosen to purchase a high-quality ADE brand product, which combines intelligent functions with an exceptional design. This set of body analyser scales enables you to measure and get to know yourself and your body better.

The wealth of experience of the ADE brand ensures a high technical standard and proven quality.

We wish you much enjoyment while weighing yourself and discovering your values!

Your ADE Team

General information

About this operating manual

This operating manual describes how to operate and maintain your product safely. Keep this operating manual in a safe place for future reference. If you pass this product on to someone else, pass on this operating manual as well.

Failure to observe this operating manual may result in injuries or damage to your product.

Explanation of symbols

This symbol combined with the word DANGER warns of serious injuries.

This symbol combined with the word WARNING warns of moderate and mild injuries.

This symbol combined with the word NOTE warns of material damage.

This symbol indicates additional information and general notes.

Table of Contents

General information	18
Scope of delivery	19
Safety	20
Intended use	21
Overview	22
Start-up	23
This is how your body analyser scales work.....	24
Measuring your weight	24
Entering Personal Data	25
Performing Body Analysis	25
Automatic user recognition	30
Cleaning	30
Fault/Remedy.....	31
Technical data	31
Declaration of Conformity.....	32
Warranty	32
Disposal.....	33

Scope of delivery

- Bathroom scale
- 3x AAA batteries, 1.5 V 
- Operating Manual

Safety

- This product may be used by children aged 8 and above, as well as by persons with reduced physical, sensory, or mental capacities or lacking experience and knowledge if they are supervised or have been instructed in the safe use of the product and understand the risks that may arise.
- Children must not play with the product.
- Cleaning and user maintenance must not be carried out by children, unless they are 8 years of age or older and are being supervised.
- Batteries, when swallowed, can be life-threatening. You should therefore keep the batteries and scale out of reach of small children. If a battery has been swallowed, medical assistance must be sought immediately.

DANGER for children

- Keep packaging material away from children. There is a risk of choking if swallowed. Also pay attention that children do not pull the packaging material over their heads.

RISK of fire/burning and/or explosion

- There is a risk of explosion if the battery is replaced improperly. Replace the battery using only the same or an equivalent type of battery. Pay attention to the "Technical data".
- Batteries must not be recharged, reactivated by other means, dismantled, thrown into fire, or short-circuited.
- Avoid contact with skin, eyes, and mucosa. If contact is made with the battery fluid, rinse the affected areas immediately with plenty of clean water and see a doctor straight away.

WARNING of risk of injury

- Do not step onto the outer edge of any one side of the scale. Risk of tipping!
- Do not step onto the scale with wet feet or when its surface is damp. Risk of slipping!
- The weighing surface is made of glass. If you fall down, there is a risk of the glass breaking and thus a risk of injury from glass splinters.

BEWARE of material damage

- Place the scale on a stable, even surface.
- Do not overload the scale. Pay attention to the “Technical data”. Place yourself on the scale exclusively for weighing purposes.
- Protect the scale from strong shocks and vibrations.
- Remove the batteries from the scale when they are depleted or when you are no longer using the scale. You will thus avoid damage that may result from leakage of battery acid.
- When inserting and changing batteries, pay attention to the correct polarity (+/-).
- Do not subject the batteries to any extreme conditions, by – for example – storing them on radiators or under direct sunlight. Increased risk of leakage!
- Clean the battery and device contacts, as needed, before inserting the battery.
- Do not make any modifications to the product. Only allow repairs to be carried out by a specialist workshop. Repairs carried out improperly can considerably endanger the user.
- There are feet on the underside of the scale. Since floor surfaces are made from various materials and are treated with diverse care products, it cannot be completely excluded that some of these substances contain ingredients that may adversely affect and soften the feet of the device. If necessary, place a non-slip pad under the scale.

Intended use


The bathroom scale is designed for the weighing of people in the home. It is not suitable for use in clinics, medical practices, etc. Likewise, the body analysis function is designed only for private use. The scales are not suitable for use in the professional sporting sphere. The maximum carrying capacity of 200 kg must not be exceeded.

Do not use these scales if you ...

- ... are pregnant,
- ... have a cardiac pacemaker,
- ... are a dialysis patient,
- ... suffer from oedema,
- ... have implants, e.g. artificial limbs, metal plates, etc.

If in doubt, consult your GP.

Overview

- 1 Buttons ◀ ▶
 - Set measurement unit
 - Select memory location
 - Decrease and increase data
- 2 Button ⏻
 - Enter setting mode
 - Save entries
- 3 Display
- 4 Battery compartment (bottom)
- 5 Weighing surface with integrated sensors

Start-up

For use, you need 3 batteries of type AAA, 1.5 V (included). The battery compartment is located on the underside of the device.

1. Open the battery compartment by pressing the locking tab in the direction of the arrow and removing the battery compartment cover.
2. Insert the batteries as indicated at the bottom of the battery compartment. Pay attention to the correct polarity (+/-).
3. Replace the battery compartment cover.
The locking tab must snap audibly into place.

Suggestions for changing the battery:

- A change of batteries becomes necessary when "Lo" appears on the display, or when nothing at all is displayed.
- Only use the type of battery that is specified in the "Technical data".
- Always change all batteries at the same time.
- Dispose of old batteries in an environmentally appropriate way (see "Disposal").

Choosing the unit of measurement

On delivery, the unit of measurement is set to "kg" (kilogramm) by default. However, the measurement result can also be displayed in the English unit „lb" (pound) or „st: lb" (stone pound). .

1. Step or press once briefly with little force on the weighing surface.
The display shows "0.0" and the set measurement unit, e.g. "lb".
2. Select the desired unit of measurement using button Ⓢ or Ⓣ.
Note the indication on the display.

When you enter your body data later, you must enter your height in centimeters („cm") for „kg" and for „lb" and „st: lb". These measurement units are permanently linked.

Setting up the scale

- Place the scale on an even, firm and dry surface. Do not place the scale on a carpet.

This is how your body analyser scales work

Your new body analyser scales work using so-called BIA-technology (Bioelectric Impedance Analysis). This involves an undetectable, completely harmless and weak current flowing through your body.

The measurement is performed via the sensors installed in the weighing surface.

- The current flows relatively unhindered through the fluid contained in muscle mass.
- In fat mass, the resistance is greater. This resistance is measured.

Together with your personal data (sex, weight, age, and height), the scales calculate the following information about your body composition:

- Weight in the set measurement unit
- Body fat content in % (measuring range 5-85 %)
- Body water in % (measuring range 25-75 %)
- Muscle mass in % (measuring range 25-75 %)
- Bone mass in kg (measuring range 0.5-8 kg)
- Visceral fat (measuring range 1-38)
- Basal metabolic rate BMR (calorie consumption at basal metabolic rate)
- Body mass index BMI (measure of body weight)

The scales are equipped with 8 memory cells in total.

Measurement errors

Measurement or calculation errors may arise for the following groups of people during the body analysis:

People who

... spend 10 hours or more per week playing sport,

... are professional athletes, bodybuilders or endurance athletes,

... are over 80 years old,

... are less than 10 years old.

The “normal weighing function” (see next chapter) is not affected by these considerations.

Measuring your weight

For normal weighing – without body analysis – you do not have to input any data.

- Position yourself on the scale. Stand still and do not hold on to anything. Shortly afterwards, the measured weight is displayed. A few seconds after you have left the weighing surface, the scale switches off automatically.

Entering Personal Data

- Proceed quickly when entering the data. Otherwise, the scale will switch to normal weighing mode or even switch off completely and you will have to repeat the procedure.
- If you press and hold the ⏮ or ⏭ buttons, the displayed values will advance more quickly.
- Data entered can be overwritten by re-entering it.

1. Step or press with some force on the weighing surface.
The display 2 shows "0.0" and the set measurement unit, e.g. "kg".
2. Press the  button .
3. Select the desired memory location with the ⏮ or ⏭ button.
4. Press the  button.
5. Select your gender using button ⏮ or ⏭.
6. Press .
7. Enter your height using button ⏮ or ⏭.
8. Press .
9. Enter your age using button ⏮ or ⏭.
10. Finally, press . This completes the entry of the required data.
The display shows „0" (depending on the unit of measurement), and you can start your first body analysis if you wish (see next chapter).
 - If you press  and then ⏮ or ⏭, you can switch the units of measurement.
 - If you do nothing else, the scale will switch off automatically after a few seconds.

Performing Body Analysis

- Body analysis is only possible if you have entered your data , see the previous chapter.
- Refer to the "Weighing Tips" section.


1. Step or press with some force on the weighing surface .
2. Briefly press the  button.
3. Press the ⏮ or ⏭ button to select the memory location where you have stored your personal body data.
4. Wait until the memory location number stops flashing.
5. Step barefoot onto the scale. Stand still and do not hold on to anything.
There must be no contact between the feet, legs, calves and thighs. Otherwise the measurement will be inaccurate.

Shortly afterwards, the measured weight is displayed. Afterwards, the display will show “o” and the body analysis will be performed.

Continue to stand still.

The measured values are then displayed on the display (example illustration):

1st view

2nd view

Symbol/Display	Meaning
62.40	Your body weight in kg / lb / st: lb
KCAL 1382	Your daily calorie requirement in basal metabolic rate
22.9	Your Body Mass Index (numerical value and color scale)
38 BODY age	Your determined biological age
FAT	Your body fat content in %
	Your body water in %
	Your muscle mass in %
	Your bone mass in kg
V-FAT	Your visceral fat

Note on Biological Age:

Based on the measured electrical resistance during weighing, algorithms estimate or determine your biological age using factors such as age, gender, weight, etc.

- Higher Biological Age: High body fat percentage, low muscle mass, poor metabolism.
- Lower Biological Age: Good muscle mass, low body fat percentage, active metabolism.

Weighing tips

- Weighing for the body analysis can only be carried out barefoot. In doing so, your feet should only be lightly moistened. Completely dry feet can lead to false results, since the electrical conductivity will then poor.
- During the measuring process, stand still and do not hold onto anything.
- Do not weigh yourself directly after unusual physical exertion.
- Do not weigh yourself immediately after standing up. Wait approximately fifteen minutes. In this time, the water is able to distribute through your body.
- Always attempt to weigh yourself at the same time, in order to obtain comparable values.

Notes on the ascertained values

The proportion of body fat gives the amount of fat as a percentage of the total weight of the body. A blanket recommendation for this value does not exist, since this depends significantly on sex (because of different physiques) and age.

Body fat is vital for day-to-day body function. It protects organs, cushions joints, regulates body temperature, stores vitamins, and serves the body as an energy store.

For women, a body fat proportion between 10% and 13% is essential; for men, between 2% and 5% with the exception of competitive sportsmen.

The proportion of body fat is not an unambiguous indicator of health. Excessive body weight and a too-high proportion of body fat are being linked with the appearance of many diseases of affluence such as diabetes, cardiovascular illnesses, etc.

Physical inactivity and poor diet are often triggers for these diseases. There is a clear causal connection between excess weight and physical inactivity.

A considerable part of the human body consists of water. This body water is variably distributed throughout the whole body. Lean muscle tissue contains about 75% water, blood around 83% water, body fat approx. 25% water and bones have a water proportion of about 22%.

- In men approx. 60% of body mass consists of water.
- In women, that value lies around 55% (determined by a higher proportion of body fat).

Approx. two thirds of this total body water are stored in the cells and therefore called intracellular water. The other third is extracellular water.

Fat in percent

Women				
Age	Underfat	Healthy	Overfat	Obese
10-15	5.0-16.0 %	16.1-21.0 %	21.1-29.6 %	29.7-85 %
16-30	5.0-16.4 %	16.5-22.0 %	22.1-30.4 %	30.5-85 %
31-60	5.0-17.0 %	17.1-22.4 %	22.5-31.4 %	31.5-85 %
61-80	5.0-17.6 %	17.7-23.0 %	23.1-32 %	32.1-85 %

Men				
Age	Underfat	Healthy	Overfat	Obese
10-15	5.0-13.0 %	13.1-17.4 %	17.5-25.6 %	25.7-85 %
16-30	5.0-13.6 %	13.7-18.0 %	18.1-26.4 %	26.5-85 %
31-60	5.0-14.0 %	14.1-18.6 %	18.7-27 %	27.1-85 %
61-80	5.0-14.4 %	14.5-19.0 %	19.1-28 %	28.1-85 %

Body water and muscle mass in percent

Women				
Age	Low	Healthy	High	Body Muscle
10-15	< 54.0 %	54.0-59.0 %	> 59 %	> 44 %
16-30	< 53.0 %	53.0-58.0 %	> 58 %	
31-60	< 52.0 %	52.0-57.0 %	> 57 %	
61-80	< 51.0 %	51.0-56.0 %	> 56 %	

Men				
Age	Low	Healthy	High	Body Muscle
10-15	< 57.0 %	57.0-62.0 %	> 62 %	> 44 %
16-30	< 56.0 %	56.0-61.0 %	> 61 %	
31-60	< 55.0 %	55.0-60.0 %	> 60 %	
61-80	< 54.0 %	54.0-59.0 %	> 59 %	

Bone weight (in kg)

	<i>Women</i>			<i>Men</i>		
body weight	< 50	50-75	> 75	< 75	75-99	> 99
reference bone weight	~ 5.7	~ 6.0	~ 6.3	~ 6.2	~ 6.5	~ 6.9

Daily calorie requirement

	<i>Women</i>		<i>Men</i>	
Age	weight (kg)	kilokalorie (kcal)	weight (kg)	kilokalorie (kcal)
10-17	50	1265	60	1620
18-29	55	1298	65	1560
30-49	60	1302	70	1561
50-69	60	1242	70	1505

Body mass index (BMI)

The BMI indicates the ratio of height to weight and is only an approximate standard value. According to the World Health Organization (WHO), an index of 18.6 to 25 is optimal.

Since men typically have a higher percentage of muscle mass in the total body mass than women do, the age and gender play an important role in interpreting the BMI.

BMI (kg/m²)	Body weight
< 16	Severely underfat
16 - 18,4	Underfat
18,5-24,9	Healthy
25 - 29,9	Overfat
> 30	Obese

Visceral fat

normal	slightly elevated	significantly elevated
< 9	10-14	> 15

Visceral fat is the fat around the internal organs in the abdominal cavity (e.g., liver, stomach, and intestines) and can be harmful to health.

Automatic user recognition

The scales are equipped with automatic user recognition. This function can be enabled if you have saved your personal data as described in “Entering personal data” and carried out at least one body analysis.

As soon as you step on the scales with barefoot, they switch on automatically and show you the weight result within a few seconds. When a similar weight value (difference less than 2 kg) is found in the data stored in the scales, the LCD display flashes once and shows the corresponding memory cell.

The body analysis will start shortly after.

However, automatic user recognition will not work if there is a weight difference of less than 2 kg between **two saved users**. In this case, manual memory cell selection is required.

Example

Your memory location is „P1“:

- Press .

Your memory location is „P2“:

- Press .

Cleaning

NOTE on material damage

- Do not immerse the scale in water or other liquids.
- Never use harsh, scouring or abrasive cleaning agents or the like to clean the scale. Otherwise, the surface could be scratched.
- Wipe the scale as needed with a slightly damp soft cloth.
- Disinfect the weighing surface and the sensors from time to time using a suitable disinfecting agent (available at specialist retailers and pharmacies).

Fault/Remedy

Fault	Possible cause
No function	Are the batteries empty or inserted the wrong way around?
The scale appears to be displaying wrong weighing results	- Is the scale located on a firm, even surface? - Is a device that can cause electromagnetic disturbances located in the vicinity of the scale? - Remove the batteries from the scale and wait a moment. Reinsert the batteries. The data stored up to this point will be retained.
The display shows "Err".	The maximum capacity of the scale (200 kg) has been exceeded.
"Lo" is shown on the display	The batteries are empty.
"Err2" is shown on the display	The fat content % is too low or too high and therefore cannot be displayed. Check if the entered body parameters are correct.

Technical data

Model:	BA 2500-1, BA 2500-2 body analysis scale
Battery:	3x AAA, 1.5 V 
Rated current:	60 mA
Measuring range:	max. 200 kg (min. 5 kg)
Variance:	2% +/- 100g
Indicator division	
Weight:	50 g
Body water, muscle mass, bone mass and body fat:	0.1%
Memory cells:	8
Age input:	10 – 80 years
Height:	100 – 220 cm
Scale dimensions:	approx. 320 x 280 x 25 mm

Recommended ambient conditions:

Temperatur: 0-40 °C /32-104 °F

Relative humidity: max 85 %

Distributor:
(no service address)

GRENDS GmbH
Stahltwiete 23
22761 Hamburg, Germany

Our products are constantly being developed and improved. For this reason, design and technical modifications are possible at any time.

Declaration of Conformity

GRENDS GmbH hereby declares that the body analysis scale BA 2500-x complies with EU directives 2014/30/EU and 2011/65/EU.

The full text of the EU declaration of conformity is available at the following Internet address: www.ade-germany.de/doc

This declaration will lose its validity if any modification is made to the scale without our consent.

Hamburg, March 2025
GRENDS GmbH

Warranty

GRENDS GmbH warrants for 2 years from the purchase date the cost-free rectification of defects in material and workmanship through repair or exchange.

In the event of a warranty claim, please return the product with the proof of purchase (specifying the reason for claim) to your dealer.

Disposal

Disposing of the packaging

Dispose of the packaging with similar materials. Dispose of cardboard with waste paper and of films with recyclable materials.

Disposing of the product

Dispose of the product in accordance with the regulations that apply in your country.

Devices must not be disposed of in normal household waste

At the end of its life, the product must be disposed of in an appropriate way. In this way, valuable materials contained in the device will be recycled and the burden on the environment avoided. Bring the old device to a collection point for electronic waste or to a recycling depot.

For more information, contact your local waste disposal firm or your local administrative body.

Standard and rechargeable batteries should not be disposed of along with household waste

All batteries and rechargeable batteries must be brought to a local or district collection point or to the retailer. Such batteries and rechargeable batteries can thus be disposed of in an environmentally friendly way.

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Distributeur (geen onderhoudsadres): GRENDS GmbH. Stahlwiete 23